

# Fast Synthesis Service for High-Quality DNA/RNA Oligos

## Stringent Quality Control High Guaranteed Yields Proficient Customer Support

### Fast Synthesis Service

Standard desalted oligos (genomics and 0.04  $\mu$ mol scale) ordered until 6 PM are produced overnight and are shipped the next working day (order-to-arrival time:  $\sim$ 1.5 working days!). Modified oligos undergoing a thorough purification procedure are usually produced and shipped out within 2 working days (order-to-arrival time  $\leq$ 3 working days).

### Proven Quality

As a pioneer in the field of DNA synthesis, Microsynth has gained a lot of critical know-how and experience in the synthesis of oligonucleotides since its incorporation in 1989. We are proud to supply worldwide research groups in academia, government and industry with DNA/RNA oligonucleotides meeting the highest quality standards and requirements.

### We Are Here To Help!

Do you have questions regarding your proposed application/experiment or the experimental outcome? Or do you have a special synthesis request for a tricky DNA or RNA molecule? Then please contact us and our expert staff will do its best to support you (oligo.support@microsynth.ch or 071 722 83 33). Professional customer service is of high priority for us.

### Your Benefits Include

**Large range of different synthesis scales** (up to g or even kg quantities can be realized!) – **Broad spectrum of highquality purifications** and **modifications** – Routine synthesis of DNA oligos up to 150 bases and RNA oligos up to 80 bases, or on request even

longer – **Stringent quality control** (online Trityl monitoring during synthesis on each coupling reaction and MALDI-TOF MS or analytical PAGE) – **High guaranteed yields** for the different synthesis scales/purifications – **Easy ordering** via our webshop at

www.microsynth.ch. – **Support and consulting for any special synthesis request** through our dedicated experts from the production team (no waiting on a telephone hotline!)

#### Welcome Offer

New Customers can test DNA Oligos at not charge or at reduced conditions!

Email us at [oligo.support@microsynth.ch](mailto:oligo.support@microsynth.ch)

Receive your coupon for 4 free DNA genomics oligos or a customized trial discount

#### How to order?

Login our webshop [www.microsynth.ch](http://www.microsynth.ch)

Select the blue domain „DNA/RNA Synthesis“, then click on „DNA“, „siRNA“

„single-stranded RNA“ or „2'-O-Methyl RNA“ and follow the further instructions

#### Need More Information?

Call us at +41 71 722 83 33 or

Email us at [oligo.support@microsynth.ch](mailto:oligo.support@microsynth.ch)

#### Important - Please Note!

- *Special price for large-volume orders (e.g.  $\geq$ 100 oligos or  $\geq$ 12 siRNAs) possible.*
- *Click chemistry applicable for the synthesis of more complex DNA/RNA oligos.*
- *Premium oligos can be produced for a wide range of molecular diagnostic applications.*
- *Customer-specific delivery of oligos in tubes, plates, or barcoded storage tubes possible.*

*Please contact us for a specific offer.*