

Fast Synthesis Service for High-Quality siRNAs

High Quality - all siRNAs are 100% quality controlled
Fast Delivery - shipment within 24 hours
Proficient Customer Support

High Quality

- 25 years experience in DNA synthesis and 10 years in RNA synthesis
- All siRNAs are 100% quality-controlled by MALDI-TOF MS before shipment
- Microsynth is ISO certified according to 9001:2008.

Fast Delivery

Shipment within 24 hours for single siRNA orders if your order is placed before 10 AM (short note by phone is necessary) Routine shipment within 48 hours for single siRNA orders

Guaranteed Yields

High yields are delivered for the different synthesis scales (e.g. guaranteed yield of 9 OD260 ~ 21 nmol for desalted siRNAs synthesized in the 0.04 μmol scale). Most often delivered amount significantly exceeds the minimum guaranteed yield.

Your Benefits Include

All siRNA duplexes are shipped **ready-to-use in your experiments** (fully deprotected, liquid and annealed at 40 μM) – A **free of charge siRNA design service** is available – Classical siRNA duplexes (2x21-mers) and longer duplexes up to 27-mers are available at **3 different purification grades (desalted, HPLC and PAGE)**

– Pre-synthesized **positive control siRNAs for Lamin, GFP or Luciferase** as well as **negative control siRNAs** are available on stock (4 nmol aliquotes) – **Guaranteed knock down** if you use our free siRNA design service (ask us for more detailed information) – **Easy ordering** of siRNAs via our webshop (control siRNAs via email

at oligo.support@microsynth.ch) – **Attractive prices for large-volume orders** including more than 12 siRNAs per one single order (please ask us for a specific offer) – **Support and consulting for any special synthesis request** through our dedicated expert staff from the production

Welcome Offer

- Test the siRNA synthesis service at no charge!
- Email us at info@microsynth.ch
- Receive 1 free siRNA (per research group)

How to order?

- Login our webshop www.microsynth.ch
- Select the blue domain „DNA/RNA Synthesis“
- Click on “siRNA” and follow the further instructions

Need More Information?

- Call us at +41 71 722 83 33 or
- Email us at info@microsynth.ch

Important - Please Note!

- Large-volume orders (> 12 siRNAs) can be delivered in 96-well plates fully deprotected, liquid and annealed (40 μM).
- If you need modified (5' or 3') siRNA or RNA synthesis, please contact us. Microsynth offers a broad spectrum of modifications (incl. click chemistry) and purification grades. RNA oligos can be routinely produced up to 80 RNA bases or on request even longer.